

Dry Clinical Chemistry Analysis

# SPOTCHEM D-Concept™

SD-9810,3810,4810,4820



A test system leading the way with  
excellent scalability and stylish design

Dry Clinical Chemistry Analysis

# SPOTCHEM D-Concept™

SD-9810,3810,4810,4820



A test system  
with excellent scalability  
and stylish design

## **Proprietary technology molded into a compact form**

ARKRAY's dry chemistry/immunological  
analysis technology, condensed.

## **The flexibility to combine the units you want**

Separate operation and measurement  
units: measurement units can be added to  
fit your needs.

## **An innovation in stylish design**

A novel, smart and functional design for  
clinical laboratory equipment



Calibration by magnetic card

Operation unit

## SPOTCHEM™ D-00 SD-9810

Once connected, measurement units can be collectively operated through this, the SPOTCHEM™ D-00. It can readily be calibrated just by loading a magnetic card into SPOTCHEM™ D-00. Simple card-swipe calibration.



Just set the Reagent pack, and you are ready to test.



Compact reagent pack

Mobile immuno-nephelometry analysis Measurement unit

## SPOTCHEM™ D-01 SD-3810

Accurate measurement of CRP, rheumatoid factor, urinary-microalbumin and urine creatinine. Disposable reagent packs remove the need for waste liquid treatment meaning extra safety for you.



Electrolyte measuring reagent



Biochemistry measuring reagent

Dry chemistry clinical analysis measurement unit

## SPOTCHEM™ D-02 SD-4810

Choose from 22 biochemistry items and 3 electrolyte items. Simultaneous measurement of up to 12 biochemistry items and 3 electrolyte items. Set the reagent, push the start key and you are ready to test. The solid phase reagent is individually-packaged, so you can perform measurements at your own pace.



Choose necessary items to measure



Built-in centrifuge



Choose necessary items to measure

Dry Clinical Chemistry analysis Measurement unit

## SPOTCHEM™ D-03 SD-4820

Choose from 22 biochemistry items. Simultaneous measurement of up to 12 biochemistry items. A built-in centrifuge allows for whole-blood sample testing without pretreatment (except for Hb).

SPOTCHEM D-00 (SD-9810) specifications	
Number of connectable units	3(Max.)
Warm up time	Approx.1 minute
Screen	Full- color LCD 320×240 dpi
Printer	58 mm thermal printer paper
Data storage capacity	100 measurement results for each connected unit
External output	RS-232C×1 port
Measurement environment	Temperature:10-30℃,Humidity:20-80% R.H.(with no condensation)
Dimensions	408(W)×330(D)×103(H)mm
Weight	Approx. 5 kg
Power supply voltage (main device)	AC100-240V, 50/60Hz
Power input	Up to 300VA

SPOTCHEM D-01 (SD-3810) specification	
Sample types	whole blood, serum, plasma, urine (sample types vary depending on reagent)
Reagent	SPOTCHEM D reagent pack
Measurement items	CRP, RF, $\mu$ Alb/Cre(Urine)
Measurement principle	Latex immunological nephelometry (LPIA method)、transmitted light volume by coloring substance
Measurement wavelengths	465nm, 550nm, 660nm
Loaded reagents	2
Processing speed	9-12 min. per sample, 14-18 min. required for simultaneous dual sample measurement
Minimum sample volume	50 $\mu$ L
Sample container	All-In-One Sample container
Warm up time	Approx 5 min, (From a room temperature of 25℃)
Measurement environment	Temperature:10-30° C,Humidity:20-80% R.H.(with no condensation)
Outer dimensions	408(W)×330 (D)×132(H) mm
Weight	Approx. 8kg
Power supply voltage	DC24V, 2A (supplied from operating unit)

SPOTCHEM D02 (SD-4810) specification	
Sample types	Biochemistry:serum, plasma (whole blood for Hb) Electrolyte:whole blood, serum, plasma and urine
Reagent	SPOTCHEM D single reagent, SPOTCHEM D multi reagent, SPOTCHEM D electrolyte plate
Measurement items	Biochemistry:ALB,ALP,ALT (GPT),AMY,AST (GOT),Ca, CK,CRE,FRA,GGT,GLU,Hb,HDL,IP,LD,Mg,T-BIL,TC,TG, TP,UA,UN (22 items in total)Electrolyte: Na,K,Cl Simultaneous measurement for a maximum of 15 items. (Biochemistry single 6 items + Biochemistry multi 6 items + Electrolyte 3 items)
Measurement principle	Biochemistry: Endpoint method and rate method by dual-wavelength reflectance Electrolyte: potentiometry method by ion-selective electrode technique
Measurement wavelength	405, 550, 575, 610, 820nm
Loadable reagent	6 single test strips, 1multi test strip and 1 electrolyte plate
Processing speed	Time required for 6 single test strips + 1 multi test strip is approx.18 min.Measurement time for Electrolyte plate is approx.4 mins
Minimum sample volume	30 $\mu$ L plus test consumption volume per measurement (Biochemistry: about 6 $\mu$ L per item Electrolyte: about 22 $\mu$ L)
Sample container	ARKRAY sample tube
Warm up time	Approx 8 min, (From a room temperature of 25° C)
Measurement environment	Temperature:10-30° C,Humidity:20-80%R.H.(with no condensation)
Outer dimensions	408 (w)×330(d)×132(h)mm
Weight	approx. 10 kg
Power supply voltage (main body)	DC24V, 2A (supplied via operation unit)

SPOTCHEM D-03 (SD-4820) specifications	
Whole blood, serum, plasma	
SPOTCHEM D single reagent, SPOTCHEM D multi reagent	
Biochemistry: ALB,ALP,ALT (GPT),AMY,AST (GOT),Ca,CK,CRE, FRA,GGT, GLU,Hb,HDL,IP,LD,Mg,T-BIL,TC,TG,TP, UA,UN (22 items in total) Simultaneous measurement with a maximum of 12 items (Biochemistry single: 6 items + Biochemistry multi: 6 items)	
Endpoint method and rate method by dual-wavelength reflectance	
6 single test strips, 1multi test strip	
Time required for 6 single test strips + 1 multi test strip is approx. 18 mins.Centrifuge time: 2.5 min.	
Cuvette(serum, plasma):30 $\mu$ L and test consumption volume per measurement (approx. 6 $\mu$ L per item) Centrifuge cup (whole blood):approx 530 $\mu$ L	
ARKRAY sample tube, ARKRAY centrifuge cup	
Approx 11 kg	
DC24V,2A (supplied from operating unit)	

\*Design and specifications may be subject to change without notice.

**arkray, inc.**

Kyoto Miyuki Bldg. 10F, 689 Takanna-cho,  
Nakagyo-ku,Kyoto 604-8153, JAPAN

TEL +81-75-662-8967 FAX +81-75-662-8973

