

Advanced Blood Glucose Test Meter

# GLUCOCARD™ X-mini plus

GT-1960



More advanced SMBG meter  
which successfully benefits  
the user with no coding,  
top specifications,  
and stylish design.

Advanced Blood Glucose Test Meter

# GLUCOCARD™ X-mini plus

GT-1960



## Design

Impressive & Stylish design

## Operation

Comprehensive display  
New user friendly functions

## Small Sample

0.3uL micro sampling

## Maximum flexibility

Auto coding  
Short fill detection  
Strip discard button  
Test strip port light  
Backlight LCD  
External output



Test strip for GLUCOCARD X-mini plus

Specification			
Test	Glucose in whole blood	Memory	500 test results
Sample Volume	Approx. 0.3uL	Operating Environment	Temperature range 10 to 40°C (50 to 104°F) Humidity 20 to 80% RH (Without condensation)
Test Strip	GLUCOCARD X-SENSOR	Dimensions	102x33x18.5mm
Measuring Range	10 TO 600mg/dL (0.6 to 33.3 mmol/L)	Weight	Approx. 30g
Measuring Time	5 seconds	※Design and specifications may be changed without prior notice	
Power source	2 X 3-volt Lithium battery (CR2032)		

arkray global business, inc.

Kyoto Miyuki Bldg. 10F, 689 Takanna-cho,  
Nakagyo-ku, Kyoto 604-8153, JAPAN  
TEL +81-75-662-8967 FAX +81-75-662-8973

